



Medline Industries Ltd
 3rd floor, Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington
 WA4 6HL
 United Kingdom
 TEL - +44 844 334 5237
 FAX - +44 844 334 5238

www.medline.com/uk



Technical Sheet

Flex One™

Reference(s)

Code	Description	Qty/Case
FO200	Disposable soft liner 2 L, white lid, without filter	66
FO201	Disposable soft liner 2 L, green lid, with filter	54
FO201PG	Disposable soft liner 2 L, green lid, with filter, pre-gelified	54
CFO200	Reusable rigid canister 2 L, white	15
CFO201	Reusable rigid canister 2 L, green	15

Description

Disposable soft liners:
 - With white lid, without filter: One-way valve on patient port
 - With green lid, with filter: Hydroscopic filter on vacuum port, one-way valve on patient port
 - With green lid, with filter, pre-gelified: Hydroscopic filter on vacuum port, one-way valve on patient port and gelling agent packed in soluble film

Reusable rigid canisters:
 - White : Graduated canister, must be used with disposable soft liner with white lid
 - Green : Graduated canister, must be used with disposable soft liner with green lid

Main materials

Disposable soft liners:
 - Liner: polyethylene
 - Lid: polypropylene
 - Gelling agent (Reference : FO201PG): Sodium polyacrylate in soluble film

Reusable rigid canisters:
 - Polycarbonate

Sterilisation

Non sterile

Pack factor/Packaging

FO200: 2 polyethylene bags with 33 ea
 FO201 & FO201PG: 2 polyethylene bags with 27 ea

Manufacturer

Medline International France SAS
 EN ISO 13485

Regulation

In compliance with Medical Device Regulation (EU) 2017/745
 MD Class: I Non sterile
 • Please refer to complete labeling and instruction for use before using the product.



© 2020 Medline Industries Inc. or one of its subsidiaries. All rights reserved.
 TDS_FLEX-ONE_EN03