



Medline Industries Ltd
 No 5 Booths Park - Chelford Road
 Knutsford - Cheshire
 WA16 8GS – England
 Phone: 0844 334 5237
 Fax : 0844 334 5238
www.medline.com.uk



Technical sheet

FOUR CANISTER TROLLEY

Reference

ORTROL4C

Description

Four canister roll stand intended to be used throughout the hospital to hold up to 4 canisters. It allows flexibility in the location and configuration of Medline suction system:

MED-SOFT, MED-RIGID, MED SEMI-RIGID and FLEX-ONE

This 5 wheels stand has been specifically designed to ensure total stability and optimal functionality of MEDLINE suction systems, which must be fixed in a vertical position during use.

The trolley is equipped with 2 wheel brakes for extra safety.

Cleaning: Reusable devices can be wiped with disinfectant.

Main materials

Nylon
 PVC
 Latex free

Sterilisation

Non sterile

Pack factor

2 units per case

Manufacturer

Prosimed Srl - EN ISO 13485
 Distributed by Medline International France in Europe

Regulation

In compliance with Medical Device Directive (EU) 93/42/CEE
 MD Class: I Non sterile
 Please refer to complete labeling and instruction for use before using the product.

