



www.medline.eu/uk



Technical Data Sheet

ACCUtouch®

Non Sterile Powder Free Synthetic Exam Glove

Reference & Size

MDS192073	MDS192074	MDS192075	MDS192076	MDS192077
XS	S	M	L	XL

Primary Material

Latex Free
Powder Free in accordance with EN455-3
Polyvinyl Chloride

Color

Clear

Surface

Smooth

Cuff Design

Beaded

Chemical additive (accelerators)

None

Powder content (per EN455-3)

Powder Content:0.12 milligrams per glove

Thickness (based on average)

Finger Tip	0.10 mm
Palm	0.08 mm
Cuff	0.06 mm

Packaging

100 gloves per box , 10 boxes per case (XL – 90 gloves per box)
Box dimensions : 240mm L x 120mm W x 60mm H
Case dimensions : 311mm L x 251mm W x 245mm H

Features:

- Flat Layer Packaging
- Color coded sizing
- Microstop Film



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-MDS19207 X
-EN03

01.INT.03.30

Length & Width Per EN455-2

MDS192073	MDS192074	MDS192075	MDS192076	MDS192077
250 mm	252 mm	250 mm	255 mm	255 mm
77 mm	86 mm	93 mm	105 mm	115 mm

Force @ Break Before Challenge Testing (per EN455-2 ≥ 3.6 N)

4.7 N

Force @ Break After Challenge Testing (per EN455-2 ≥ 3.6 N)

3.8 N

Freedom from Holes (per EN 455 AQL 1.5)

1.5 AQL

Viral Penetration

Tested and passed, in accordance with EN 374-5

Chemical Resistance

The resistance to chemicals has been assessed in accordance with EN 16523-1 Results and recommendations for use with chemicals can be obtained on request

Expiration Date

35 Months from Date of Manufacture
Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

Regulations and Quality Standards

Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

PPE Certification

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5

Storage Recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone..

Country of Origin

China

Legal Manufacturer

Medline International France SAS
5 rue Charles Lindbergh, 44110 Châteaubriant, France

