



www.medline.eu/uk



Technical Data Sheet

Medline Sterile Procedure Nitrile

Sterile Nitrile – Powder Free Glove

Reference & Size

MDS2294	MDS2295	MDS2296	MDS2297
S	M	L	XL

Primary Material

Latex Free
Sterile Powder-Free Nitrile

Color

Blue

Surface

Tip Textured

Cuff Design

Beaded cuff

Chemical additive (accelerators)

Carbamate

Powder Content (per EN455-3)

Powder Content: 0.12 milligrams per glove

Thickness (based on average)

Finger Tip	0.14 mm
Palm	0.11 mm
Cuff	0.09 mm

Packaging

50 gloves per box , 6 boxes per case (XL –50 gloves per box)
Box: 300mm L x 150mm W x 135mm H
Case: 470mm L x 285mm W x 310mm H



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-MDS229X-
EN02

01.INT.03.30

Length & Width Per EN455-2

MDS2294	MDS2295	MDS2296	MDS2297
302 mm	299 mm	301 mm	307 mm
86 mm	96 mm	105 mm	114 mm

Force @ Break Before Challenge Testing (per EN455-2)

7.8 N

Force @ Break After Challenge Testing (per EN455-2)

8.3 N

Freedom from Holes (per EN 455-1)

Pass 1.5 AQL (standard requests compliance with AQL1.5)

Viral Penetration

Tested and passed, in accordance with EN374-5

Chemical Resistance

The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978.

Results and recommendations for use with chemicals and chemo drugs can be obtained on request

Expiration Date

35 Months from Date of Manufacture
Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

Regulations and Quality Standards

"Class 1 sterile Medical Device per Directive 93/42/EEC. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

PPE Certification

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5

Storage Recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.

Country of Origin

Malaysia

Legal Manufacturer

Medline International France SAS
5 rue Charles Lindbergh, 44110 Châteaubriant, France