



www.medline.eu/uk



Technical Data Sheet

SensiCare® Free

Non sterile Nitrile – Powder Free - Accelerator Free Glove

Reference & Size

MDSFREE-XS	MDSFREE-S	MDSFREE-M	MDSFREE-L	MDSFREE-XL
XS	S	M	L	XL

Primary Material

Latex Free
Powder-Free Nitrile
Accelerators Free

Color

Blue

Surface

Finger Tip Textured

Cuff Design

Beaded cuff

Powder content (per EN455-3)

Powder Content:0.10 milligrams per glove

Thickness (based on average)

Finger Tip	0.09 mm
Palm	0.07 mm
Cuff	0.05 mm

200 gloves per box , 10 boxes per case (XL – 180 gloves per box)
Box dimensions : 240mm W x 120mm L x 77mm H
Case dimensions : 400mm W x 248mm L x 255mm H

Packaging

Features:

- MicroStop Film
- Flat Layer Packaging
- Color coded sizing



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-MDSFREE-X
-EN02

01.INT.03.30

Length & Width Per EN455-2

MDSFREE-XS	MDSFREE-S	MDSFREE-M	MDSFREE-L	MDSFREE-XL
77 mm	87 mm	96 mm	108 mm	117 mm
240 mm	251 mm	250 mm	250 mm	250 mm

Force @ Break Before Challenge Testing (per EN455-2 ≥ 6 N)

7.1 N

Force @ Break After Challenge Testing (per EN455-2 ≥ 6 N,

7.0 N

Freedom from Holes (per EN 455-1)

1.5 AQL

Viral Penetration

Tested and passed, in accordance with EN374-5

Chemical Resistance

The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978.

Results and recommendations for use with chemicals and chemo drugs can be obtained on request

Expiration Date

35 Months from Date of Manufacture

Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

Regulations and Quality Standards

Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

PPE Certification

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5

Food Contact

Tested according to European Regulation (EU) No 10/2011

Storage Recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.

Country of Origin

Malaysia

Legal Manufacturer

Medline International France SAS

5 rue Charles Lindbergh, 44110 Châteaubriant, France

