



www.medline.eu/uk



## Technical Data Sheet

# Restore<sup>®</sup> Sense

### Nitrile Powder Free Glove with Colloidal Oatmeal

#### Reference & Size

OATN32XS	OATN32S	OATN32M	OATN32L	OATN32XL
XS	S	M	L	XL

#### Primary Material

Latex Free  
 Powder-Free Nitrile  
 Contains Colloidal Oatmeal – Skin moisturiser

#### Color

Green Lemon

#### Surface

Finger Tip Textured

#### Cuff Design

Beaded

#### Chemical additive (accelerators)

Carbamates

#### Powder content (per EN455-3 ≤2mg/glove)

Powder Content: 0,52 mg/glove

#### Thickness (based in minimum)

Finger Tip	0.05 mm
Palm	0.05 mm
Cuff	0.04 mm

#### Packaging

250 gloves per box , 10 boxes per case (XL – 230 gloves per box)  
 Box dimensions : 128mm W x 240mm L x 83mm H  
 Case dimensions : 262mm W x 397mm L x 253mm H

#### Features:

- Micro-Stop Film
- Flat Layer Packaging
- Color coded sizing



**Dimensions :**  
(Length & Width per EN455-2)

OATN32XS	OATN32S	OATN32M	OATN32L	OATN32XL
76 mm	86 mm	97 mm	108 mm	115 mm
246 mm	250 mm	250 mm	250 mm	250 mm

**Force @ Break Before Challenge Testing (per EN455-2 ≥ 6 N)**

6.8 N

**Force @ Break After Challenge Testing (per EN455-2 ≥ 6 N, 13 days 70°C in an oven)**

6.2 N

**Freedom from Holes (per EN455-1)**

Pass 1.0 AQL (standard requests compliance with AQL1.5)

**Viral Penetration**

Tested and passed, in accordance with EN 374-5

**Chemical Resistance**

The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978.  
Results and recommendations for use with chemicals and chemo drugs can be obtained on request

**Expiration Date**

35 Months from Date of Manufacture  
Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

**Regulations and Quality Standards**

Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

**PPE Certification**

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5 and ISO 16604.

**Food Contact**

Tested according to European Regulation (EU) No 10/2011

**Storage Recommendations**

Protect from freezing. Avoid excessive heat. Keep dry. Product should be kept away from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.

**Country of Origin**

Malaysia

**Legal Manufacturer**

Medline International France SAS  
5 rue Charles Lindbergh, 44110 Châteaubriant, France

