



www.medline.eu/uk



Technical Data Sheet

Restore[®]

Nitrile Powder Free Glove with Colloidal Oatmeal

Reference & Size

OAT-XS	OAT-S	OAT-M	OAT-L	OAT-XL
XS	S	M	L	XL

Primary Material

Latex Free
Powder-Free Nitrile
Contains Colloidal Oatmeal – Skin moisturiser

Color

Natural White

Surface

Finger Tip Textured

Cuff Design

Beaded

Chemical additive (accelerators)

Carbamates

Powder content (per EN455-3 ≤2mg/glove)

Powder Content: 1,56 mg/glove

Thickness (based on average)

Finger Tip	0.09 mm
Palm	0.07 mm
Cuff	0.05 mm

Packaging

200 gloves per box , 10 boxes per case (XL – 180 gloves per box)
Box dimensions : 245 mm L x 125mm W x 74mm H
Case dimensions : 382mm L x 262mm W x 257mm H

Features:

- MicroStop Film
- Flat Layer Packaging
- Color coded sizing



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-OAT-X-EN02

01.INT.03.30

Dimensions :
(Length & Width per EN455-2)

OAT-XS	OAT-S	OAT-M	OAT-L	OAT-XL
78 mm	87 mm	99 mm	108 mm	117 mm
246 mm	250 mm	242 mm	252 mm	255 mm

Force @ Break Before Challenge Testing (per EN455-2 \geq 6 N)

6.3 N

Force @ Break After Challenge Testing (per EN455-2 \geq 6 N, 13 days 70°C in an oven)

6.1 N

Freedom from Holes (per EN455-1)

Pass 1.0 AQL (standard requests compliance with AQL1.5)

Viral Penetration

Tested and passed, in accordance with EN 374-5

Chemical Resistance

The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978.

Results and recommendations for use with chemicals and chemo drugs can be obtained on request

Expiration Date

35 Months from Date of Manufacture
Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

Regulations and Quality Standards

Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

PPE Certification

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5 and ISO 16604.

Food Contact

Tested according to European Regulation (EU) No 10/2011

Storage Recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be kept away from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.

Country of Origin

Malaysia

Legal Manufacturer

Medline International France SAS
5 rue Charles Lindbergh, 44110 Châteaubriant, France

