

Standard nitrile



MediGuard Special Nitrile

Item number	MGSNXS MGSNS MGSNM MGSNL MGSNXL
Sizes	XS – XL
Colour	Violet blue
Texture	Fingertips
Force at break (Unchallenged)	7.0 N
AQL	1.0
Cuff length (Size M)	240 mm
Thickness (Fingers; Palm; Cuff)	0.010; 0.08; 0.07 (± 0.03) mm
Packaging	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



SensiCare® Silk

Item number	MDS9083 MDS9084 MDS9085 MDS9086 MDS9087
Sizes	XS – XL
Colour	Violet blue
Texture	Fingertips
Force at break (Unchallenged)	6.7 N
AQL	1.0
Cuff length (Size M)	244 mm
Thickness (Fingers; Palm; Cuff)	0.09; 0.07; 0.05 mm
Packaging	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



SensiCare® Ice

Item number	486800 486801 486802 486803 486804
Sizes	XS – XL
Colour	Violet blue
Texture	Fingertips
Force at break (Unchallenged)	8.1 N
AQL	1.0
Cuff length (Size M)	249 mm
Thickness (Fingers; Palm; Cuff)	0.09; 0.07; 0.05 mm
Packaging	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



Exam Gloves Ordering Guide

Special nitrile



SensiCare® Free Accelerator-Free

Specialty	Accelerator-Free
Item number	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL
Sizes	XS – XL
Colour	Blue
Texture	Fingertips
Force at break (Unchallenged)	8.2 N
AQL	1.0
Cuff length (Size M)	243 mm
Thickness (Fingers; Palm; Cuff)	0.09; 0.07; 0.05 mm
Packaging	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



Restore® Sense Colloidal Oatmeal Inner Coating

Item number	OATN32XS OATN32S OATN32M OATN32L OATN32XL
Sizes	XS – XL
Colour	Lemon green
Texture	Fingertips
Force at break (Unchallenged)	6.8 N
AQL	1.0
Cuff length (Size M)	244 mm
Thickness (Fingers; Palm; Cuff)	0.05; 0.05; 0.045 mm
Packaging	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs



Generation Pink Nitrile Breast cancer awareness

Item number	PINKE2500 PINKE2501 PINKE2502 PINKE2503 PINKE2504
Sizes	XS – XL
Colour	Pink
Texture	Fingertips
Force at break (Unchallenged)	6.7 N
AQL	1.0
Cuff length (Size M)	240 mm
Thickness (Fingers; Palm; Cuff)	0.09; 0.07; 0.05 mm
Packaging	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs



Sterile Procedure Nitrile Sterility

Item number	MDS2294 MDS2295 MDS2296 MDS2297
Sizes	S – XL
Colour	Blue
Texture	Fingertips
Force at break (Unchallenged)	7.8 N
AQL	1.0
Cuff length (Size M)	302 mm
Thickness (Fingers; Palm; Cuff)	0.14; 0.011; 0.09 mm
Packaging	S – L: 50 pairs/bx, 6 bx/cs XL: 50 pairs/bx, 6 bx/cs

Long-cuff nitrile



SensiCare® Extended Cuff

Item number	MDSEXTXS MDSEXTS MDSEXTM MDSEXTL MDSEXTXL
Sizes	XS – XL
Colour	Blue
Texture	Fingertips
Force at break (Unchallenged)	12.4 N
AQL	1.0
Cuff length (Size M)	280 mm
Thickness (Fingers; Palm; Cuff)	0.14; 0.10; 0.08 mm
Packaging	XS – L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs



SensiCare® Pro

Item number	PRO400XS PRO400S PRO400M PRO400L PRO400XL
Sizes	XS – XL
Colour	Blue
Texture	Fully textured
Force at break (Unchallenged)	20.4 N
AQL	1.0
Cuff length (Size M)	412 mm
Thickness (Fingers; Palm; Cuff)	0.23; 0.18; 0.11 mm
Packaging	XS – L: 50 ea/bx, 10 bx/cs XL: 50 ea/bx, 10 bx/cs

Vinyl



Eco Vinyl

MGE323
MGE324
MGE325
MGE326
MGE327

Item number

Sizes

XS - XL

Colour

Clear

Texture

Smooth

Force at break (Unchallenged)

4.4 N

AQL

1.5

Cuff length (Size M)

245 mm

Thickness

(Fingers; Palm; Cuff)

0.11; 0.09; 0.07 (± 0.03) mm

Packaging

XS - L: 150 ea/bx, 10 bx/cs
XL: 130 ea/bx, 10 bx/cs



ACCUtouch®

MDS192074
MDS192075
MDS192076
MDS192077

S - XL

Clear

Smooth

4.7 N

1.5

250 mm

0.10; 0.08; 0.06 mm

S - XL: 100 ea/bx, 10 bx/cs



Curad® 3G

CUR8233
CUR8234
CUR8235
CUR8236
CUR8237

XS - XL

White

Smooth

3.7 N

1.5

248 mm

0.10; 0.08; 0.06 mm

XS - L: 100 ea/bx, 10 bx/cs
XL: 90 ea/bx, 10 bx/cs

Latex



Curad® Latex

CURE8103
CURE8104
CURE8105
CURE8106
CURE8107

XS - XL

White

Smooth

7.5 N

1.5

249 mm

0.13; 0.11; 0.09 mm

XS - L: 100 ea/bx, 10 bx/cs
XL: 90 ea/bx, 10 bx/cs

We provide
you the right
gloves for the
right task.

Wire box holders



Single box holder

MDS191094

Item number

Material

Coated wire

Colour

White

Packaging

10 holders per case



Double box holder

MDS192094B

Coated wire

White

6 holders per case



Triple box holder

MDS193094

Coated wire

White

10 holders per case

Plexiglass box holders



Single box holder

INTL191096

Item number

Material

Plexiglass

Colour

Clear

Packaging

12 holders per case



Triple box holder

INTL193096

Plexiglass

Clear

6 holders per case