



# EXAM GLOVES GUIDE

Product Name	Item Number	Size	Color	Texture	Thickness			Cuff Length	Description	Packaging
					Fingers	Palm	Cuff			
 <b>Restore Nitrile**</b> 	OAT-XS OAT-S OAT-M OAT-L OAT-XL	X-Small Small Medium Large X-Large	White	Textured Finger Tips	0.09 mm	0.07 mm	0.05 mm	241 mm	Nitrile, Non-Powdered, Non-Sterile, With colloidal oatmeal inner coating	<b>XS - L:</b> 200 ea/Box, 10 Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>SensiCare® Free**</b> 	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL	X-Small Small Medium Large X-Large	Dark Blue	Textured Finger Tips	0.10 mm	0.10 mm	0.08 mm	244 mm	Nitrile, Non-Powdered, Non-Sterile, Accelerator-Free	<b>XS - L:</b> 200 ea/Box, 10 Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>Generation Pink Nitrile**</b>	PINKE6800 PINKE6801 PINKE6802 PINKE6803 PINKE6804	X-Small Small Medium Large X-Large	Pink	Textured Finger Tips	0.09 mm	0.07 mm	0.05 mm	244 mm	Nitrile, Non-Powdered, Non-Sterile	<b>XS - L:</b> 200 ea/Box, 10/Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>SensiCare® Ice**</b> 	486800 486801 486802 486803 486804	X-Small Small Medium Large X-Large	Blue	Textured Finger Tips	0.09 mm	0.07 mm	0.05 mm	241 mm	Nitrile, Non-Powdered, Non-Sterile	<b>XS - L:</b> 200 ea/Box, 10 Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>SensiCare® Silk**</b> 	MDS9083 MDS9084 MDS9085 MDS9086 MDS9087	X-Small Small Medium Large X-Large	Dark Blue	Textured Finger Tips	0.09 mm	0.06 mm	0.05 mm	245 mm	Nitrile, Non-Powdered, Non-Sterile	<b>XS - L:</b> 200 ea/Box, 10 Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>MediGuard Blue**</b> 	EM5050 EM5051 EM5052 EM5053 EM5054	X-Small Small Medium Large X-Large	Light Blue	Textured Finger Tips	0.13 mm	0.07 mm	0.05 mm	246 mm	Nitrile, Non-Powdered, Non-Sterile	<b>XS - L:</b> 200 ea/Box, 10 Boxes/Case <b>XL:</b> 180 ea/Box, 10 Boxes/Case
 <b>SensiCare® Pro**</b> 	PRO400XS PRO400S PRO400M PRO400L PRO400XL	X-Small Small Medium Large X-Large	Blue	Fully Textured	0.23 mm	0.18 mm	0.11 mm	405 mm	Nitrile, Non-Powdered, Non-Sterile	<b>XS - XL:</b> 50 ea/Box, 10 Boxes/Case
 <b>Aloetouch® Long Cuff**</b>	MDS195184 MDS195185 MDS195186 MDS195187	Small Medium Large X-Large	Green	Fully Textured	0.19 mm	0.17 mm	0.09 mm	297 mm	Nitrile, Non-Powdered, Non-Sterile, Aloe Coating Inside	<b>S - L:</b> 50 ea/Box, 10 Boxes/Case <b>XL:</b> 50 ea/Box, 10 Boxes/Case
 <b>SensiCare® Long Cuff**</b> 	MDS1284 MDS1285 MDS1286 MDS1287	Small Medium Large X-Large	Blue	Textured Finger Tips	0.14 mm	0.11 mm	0.05 mm	292 mm	Nitrile, Non-Powdered, Non-Sterile	<b>S - L:</b> 50 ea/Box, 10 Boxes/Case <b>XL:</b> 45 ea/Box, 10 Boxes/Case

\*\* These gloves have been tested and approved for use with chemotherapy drugs

# EXAM GLOVES GUIDE

Product Name	Item Number	Size	Color	Texture	Thickness			Cuff Length	Description	Packaging
					Fingers	Palm	Cuff			
 <b>Curad® 3G Stretch Vinyl</b> 	CUR8233	X-Small	White	Smooth	0.09 mm	0.07 mm	0.07 mm	242 mm	Stretch Vinyl Non-Powdered Non-Sterile	<b>XS - L:</b> 100 ea/Box, 10 Boxes/Case
	CUR8234	Small								
	CUR8235	Medium								
	CUR8236	Large								
	CUR8237	X-Large								
 <b>Generation PINK 3G Vinyl</b>	PINK6073	X-Small	Pink	Smooth	0.16 mm	0.08 mm	0.07 mm	247 mm	Stretch Vinyl Non-Powdered Non-Sterile	<b>XS - L:</b> 100 ea/Box, 10 Boxes/Case
	PINK6074	Small								
	PINK6075	Medium								
	PINK6076	Large								
	PINK6077	X-Large								
 <b>ACCUtouch</b>	MDS192073	X-Small	Clear	Smooth	0.15 mm	0.08 mm	0.06 mm	254 mm	Synthetic Vinyl Non-Powdered Non-Sterile	<b>XS - L:</b> 100 ea/Box, 10 Boxes/Case
	MDS192074	Small								
	MDS192075	Medium								
	MDS192076	Large								
	MDS192077	X-Large								
 <b>Eco Vinyl</b> 	MGE323	X-Small	Clear	Smooth	0.10 mm	0.07 mm	0.05 mm	252 mm	Synthetic Vinyl Non-Powdered Non-Sterile	<b>S - L:</b> 150 ea/Box, 10 Boxes/Case
	MGE324	Small								
	MGE325	Medium								
	MGE326	Large								
	MGE327	X-Large								
 <b>Curad® Latex</b> 	CURE8103	X-Small	Beige	Textured Finger Tips	0.13 mm	0.11 mm	0.09 mm	246 mm	Latex Non-Powdered Non-Sterile	<b>XS - L:</b> 100 ea/Box, 10 Box/Case
	CURE8104	Small								
	CURE8105	Medium								
	CURE8106	Large								
	CURE8107	X-Large								
 <b>Medline Procedure Sterile Latex</b>	MDS8114	Small	Cream	Fully Textured	0.17 mm	0.15 mm	0.09 mm	243 mm	Latex Non-Powdered Sterile	<b>S - XL:</b> 50 Pairs/Box, 6 Boxes/Case
	MDS8115	Medium								
	MDS8116	Large								
	MDS8117	X-Large								
 <b>Medline Sterile Procedure Nitrile**</b> 	MDS2294	Small	Blue	Textured Finger Tips	0.12 mm	0.10 mm	0.08 mm	296 mm	Nitrile Non-Powdered Sterile	<b>S - XL:</b> 50 Pairs/Box, 6 Boxes/Case
	MDS2295	Medium								
	MDS2296	Large								
	MDS2297	X-Large								



## Medline Industries Ltd

No 5 Booths Park  
Chelford Road  
Knutsford, Cheshire WA16 8GS  
Tel: +44 844 334 5237  
Fax: +44 844 334 5238  
www.medline.com/uk  
uk-customerservice@medline.com