



PROTECT

AND PERFORM.

Longer Gloves for Extra Protection
400mm Cuff, Heavy Duty Nitrile, and Excellent Grip

EXTRA PROTECTION. BETTER GRIP.

Designed for staff that use harsh chemicals in general tasks and procedures or clean, decontaminate, and disinfect surgical instruments. Our SensiCare Pro nitrile gloves provide exceptional protection, fit, and performance.

- » **All-around protection and performance**— longer 400mm cuffs, thick glove, performs well in puncture resistance tests*, chemo-tested**, and fully textured for better grip.

Ordering Information

SensiCare® Pro 400 mm Extended Cuff Nitrile Gloves

Item No.	Size	Pkg.
PRO400XS	X-Small	50 ea/bx, 10 bx/cs
PRO400S	Small	50 ea/bx, 10 bx/cs
PRO400M	Medium	50 ea/bx, 10 bx/cs
PRO400L	Large	50 ea/bx, 10 bx/cs
PRO400XL	X-Large	50 ea/bx, 10 bx/cs



Color: Blue
Length: 400mm
Material: Nitrile
Surface: Textured
AQL: 1.5

EN374-2:2014	Performance Level	AQL	Inspection Levels
 AQL < 1.5 (level 2)	level 3	<0.65	G1
	level 2	<1.5	G1
	level 1	<4.0	S4

Additional permeation information is available upon request.

BS EN 374-3:2003 Chemical Hazards  Sodium Hydroxide (NaOH) 40% level 6 (>480 minutes)	EN374-3 Permeation levels are based on breakthrough times as follows:						
	Level	1	2	3	4	5	6
Min. breakthrough times (mins)	10	30	60	120	240	480	

Results do not necessarily reflect the actual duration of protection due to external factors.

EN388:2003	4001	Abrasion Resistance 4/4
 BS EN 388 Mechanical Hazards		Blade Cut Resistance 0/5
		Tear Resistance 0/4
		Puncture Resistance 1/4



Ideal for specialty areas, including:

-  Central Sterile
-  Endoscopy
-  Pathology
-  Environmental Services
-  Pharmacy
-  Emergency Department
-  Trauma
-  Labor and Delivery
-  Lab
-  Decontamination
-  Urology

*Test data available upon request. **Do not use with Carmustine and Thiotepea.



Medline International France SAS

5 Rue Charles Lindebergh, 44110
Chateaubriant, France
medline.com | info@medline.com