



**TRUSTED FLUID  
MANAGEMENT.**

**Suction Connecting Tubing**  
Make worry-free connections

# SECURE FLOW WITH MEDLINE SUCTION TUBING

- Strong, collapse-resistant, low memory surgical tubing.
- Non-conductive connecting tubing available in 5 mm to 7 mm inner diameter.
- Kink resistant but flexible and particularly suitable for use with high vacuum.

## STERILE SUCTION TUBING<sup>(1)</sup>

ORNEX tubing incorporates an anti-slip surface which counters any tendency for the tubes to slide off the operating table. They are structured to avoid internal bottlenecks to ensure maximum flow rate. The inside ribbed female connector provides a secure fit on a wide range of surgical instruments. Sterile double wrapped individual peel pouches.



### ORNEX STERILE

with ribbed female connectors and pre-attached male connector.

Code	ID	Length	Pack factor	NHS
ORNEX56A	5 mm	1.80 m	50/cs	FDR297
ORNEX510		3 m	30/cs	FDR296
ORNEX66A	6 mm	1.80 m	45/cs	FDR298
ORNEX610		3 m	25/cs	FDR280
ORNEX76A	7 mm	1.80 m	40/cs	FDR281
ORNEX710		3 m	20/cs	FDR279
ORNEX712		3.6 m	20/cs	FDR299



### ORNEX VC STERILE

with ribbed female connectors and pre-attached vacuum connector.

Code	Type	ID	Length	Pack factor
ORNEX56VC	F/F+MVC	5 mm	1.80 m	50/cs
ORNEX510VC	F/F+MVC		3 m	50/cs
ORNEX66VC	F/MVC	6 mm	1.80 m	50/cs
ORNEX610VC	F/MVC		3 m	50/cs



» ID means inner diameter of the tubing

» F means Female connector fits suction machines and male connections for instrumentation.

» M means Male connector push fits into the F connector and is normally used to connect to the standard funnel fitter on suction catheters. Can also be used to connect lengths of tubing together for longer lengths.



# NON-STERILE SUCTION TUBING<sup>(2)</sup>

## ORNEX WITH MED-STOP

Suction tubing with MED-STOP for removable tracheal aspiration.


Code	ID	Length	Med-Stop	Pack factor
ORNEX66VCA	6 mm	1.80 m	Adult	30/cs
ORNEX66VCP	6 mm	1.80 m	Paediatric	30/cs
ORNEX610VCA	6 mm	3 m	Adult	30/cs
OR853VCA	MED-STOP only		Adult	200/cs
OR854VCP	MED-STOP only		Paediatric	200/cs



## ORNEX ROLL TUBES

Non-sterile tubing in convenient 30 meter roll. Easy to cut to required length. Cut-to-size female connectors are located at 1.8 m intervals. Sixteen male connectors are included in the box.


Code	ID	Length	Pack factor
ORNEX5100	5 mm	30 m roll = 1.80 m x 16 suction tubes	1/cs
ORNEX6100	6 mm		1/cs
ORNEX7100	7 mm		1/cs



## NBT BUBBLE ROLL TUBES<sup>(3)</sup>

Universal suction bubble tubing. Easy to cut through the bubble to required length and/or tapered end to fit any suction catheter or suction tip.

Code	ID	Length	Pack factor	NHS
NBT5100	5 mm	30 m roll = 0.90 m x 33 suction tubes	1/cs	FDR295
NBT7100	7 mm		1/cs	FDR278



# ACCESSORIES<sup>(2)</sup>

## PLASTIC CONNECTORS

Non-sterile plastic connectors are durable and shatter-resistant. Light in weight and sufficiently transparent to allow easy observation of flow. Single use. Supplied in convenient dispenser cartons clearly labeled by type.



OR356


Code	Description	Pack factor
OR350A	<b>Straight Connector</b> Fits five sizes of tubing. Connect same size tubes or use as a reducer for any combination of 5 to 11 mm tubing.	500/cs
OR351	<b>Sims Connector</b> Fits tubing with a diameter from 3 to 6 mm in seven size increments, graduated connector adapts from 3 to 6 mm for a broad array of end user applications.	500/cs
OR352	<b>“T” connector</b> Used to join three sections of tubing. Fir tree connector for tubing connections from 5 to 11 mm.	150/cs
OR356	<b>“Y” Connectors</b> Fits six sizes of tubing from 6 mm to over 15 mm. It connects the same or different size tubes.	150/cs



OR350A



OR351



OR352

<sup>(1)</sup>Those products are class I sterile medical devices.



<sup>(2)</sup>Those products are class I non-sterile medical devices.



<sup>(3)</sup>Those products are class I non-sterile medical devices.



Before use, consult instructions for use and precautions on the corresponding labeling.



### Medline Industries Ltd

No 5 Booths Park  
Chelford Road  
Knutsford, Cheshire WA16 8GS  
England  
Tel: +44 844 334 5237  
Fax: +44 844 334 5238  
www.medline.eu/uk  
uk-customerservice@medline.com