
DATA SHEET

ESE1550EU
Disposable electrode

Pictures:



ESE1550EU Disposable electrode is fitted with high frequency surgical equipment for cutting or coagulation in surgery products; it has reasonable structure, clean appearance, excellent conductivity and easily operation. It is used in general surgery, orthopedics, cardiology, gynecology, urology, neurosurgery, dermatology, otolaryngology, abdominal surgery, plastic surgery and other surgeries for cutting and coagulation. This product is suitable for the research department of medical and health institutions, operating room and other units with high frequency surgical equipment.

TECHNICAL FEATURES	
MODEL	**ESE1550EU**
DESCRIPTION	**Disposable electrode，outer diameter Φ2.36mm**
MATERIAL	**303F + DOW 470**
LENGTH	**50.8±1mm**
ELECTRODE	**Disposable electrode**
STERILISATION	**EO sterilization**
COMPLIANCE TO ELECTRICAL SAFETY STANDARDS	**EN 60601-1 EN 60601-2-2**
COMPLIANCE TO DIRECTIVE	**Directive 93/42/EEC. Class IIB device**

PACKAGING	
PRIMARY PACKAGING	**85g Pouch**
SALE PACKAGING	**200 units in carton box 430mm(L)X155mm(W)X180mm(H)**
STORAGE	**Temperature -40 ~ +50 ℃**
EXPIRY	**3 years from production date if correctly stored**